

SNARES



Amplatz Goose Neck™
Snare Kit

Amplatz Goose Neck™

Snare Kit

Supporting successful outcomes.

Reference Number	Loop Diameter (mm)	Snare Length (mm)	Catheter Size (F)	Catheter Length (cm)
GN500	5	120	4	102
GN1000	10	120	4	102
GN1001	10	65	4	48
GN1500	15	120	6	102
GN2000	20	120	6	102
GN2500	25	120	6	102
GN2501	25	65	6	48
GN3000	30	120	6	102
GN3500	35	120	6	102

Amplatz Goose Neck™

MicroSnare Kit

Supporting successful outcomes.

Reference Number	Loop Diameter (mm)	Snare Length (mm)	Catheter Size (F)	Catheter Length (cm)
SK200	2	175	2.3-3	150
SK201	2	200	2.3-3	175
SK400	4	175	2.3-3	150
SK401	4	200	2.3-3	175
SK700	7	175	2.3-3	150
SK701	7	200	2.3-3	175