



# STENTS

EverFlex™  
Self-expanding Peripheral  
Stent System



Visi-Pro™  
Balloon-expandable Peripheral  
Stent System

# Protégé™ EverFlex™

## Self-expanding Peripheral Stent System

Flexibility sustained: Durability delivered.

CATHETER LENGTH		STENT SIZE		COMPATIBILITY		
Reference Number 80 cm	Reference Number 120 cm	Diameter (mm)	Length (mm)	Recommended Introducer Sheath (F)	Guidewire (in)	Recommended Lumen Size (mm)
PRP35-05-020-080	PRP35-05-020-120	5	20	6	0.035	3.5-4.5
PRP35-05-030-080	PRP35-05-030-120	5	30	6	0.035	3.5-4.5
PRP35-05-040-080	PRP35-05-040-120	5	40	6	0.035	3.5-4.5
PRP35-05-060-080	PRP35-05-060-120	5	60	6	0.035	3.5-4.5
PRP35-05-080-080	PRP35-05-080-120	5	80	6	0.035	3.5-4.5
PRP35-05-100-080	PRP35-05-100-120	5	100	6	0.035	3.5-4.5
PRP35-05-120-080	PRP35-05-120-120	5	120	6	0.035	3.5-4.5
PRP35-06-020-080	PRP35-06-020-120	6	20	6	0.035	4.5-5.5
PRP35-06-030-080	PRP35-06-030-120	6	30	6	0.035	4.5-5.5
PRP35-06-040-080	PRP35-06-040-120	6	40	6	0.035	4.5-5.5
PRP35-06-060-080	PRP35-06-060-120	6	60	6	0.035	4.5-5.5
PRP35-06-080-080	PRP35-06-080-120	6	80	6	0.035	4.5-5.5
PRP35-06-100-080	PRP35-06-100-120	6	100	6	0.035	4.5-5.5
PRP35-06-120-080	PRP35-06-120-120	6	120	6	0.035	4.5-5.5
PRP35-06-150-080	PRP35-06-150-120	6	150	6	0.035	4.5-5.5
N/A	PRP35DR-06-200-120	6	200	6	0.035	4.5-5.5
PRP35-07-020-080	PRP35-07-020-120	7	20	6	0.035	5.5-6.5
PRP35-07-030-080	PRP35-07-030-120	7	30	6	0.035	5.5-6.5
PRP35-07-040-080	PRP35-07-040-120	7	40	6	0.035	5.5-6.5
PRP35-07-060-080	PRP35-07-060-120	7	60	6	0.035	5.5-6.5
PRP35-07-080-080	PRP35-07-080-120	7	80	6	0.035	5.5-6.5
PRP35-07-100-080	PRP35-07-100-120	7	100	6	0.035	5.5-6.5
PRP35-07-120-080	PRP35-07-120-120	7	120	6	0.035	5.5-6.5
PRP35-07-150-080	PRP35-07-150-120	7	150	6	0.035	5.5-6.5
N/A	PRP35DR-07-200-120	7	200	6	0.035	5.5-6.5
PRP35-08-020-080	PRP35-08-020-120	8	20	6	0.035	6.5-7.5
PRP35-08-030-080	PRP35-08-030-120	8	30	6	0.035	6.5-7.5
PRP35-08-040-080	PRP35-08-040-120	8	40	6	0.035	6.5-7.5
PRP35-08-060-080	PRP35-08-060-120	8	60	6	0.035	6.5-7.5
PRP35-08-080-080	PRP35-08-080-120	8	80	6	0.035	6.5-7.5
PRP35-08-100-080	PRP35-08-100-120	8	100	6	0.035	6.5-7.5
PRP35-08-120-080	PRP35-08-120-120	8	120	6	0.035	6.5-7.5
PRP35-08-150-080	PRP35-08-150-120	8	150	6	0.035	6.5-7.5
N/A	PRP35DR-08-200-120	8	200	6	0.035	6.5-7.5

# Protégé™ GPS™

## Self-expanding Peripheral Stent System

Delivering BIGG results.

Reference Number	Unconstrained Stent Diameter (mm)	Unconstrained Stent Length (mm)	Rec. Lumen Size (mm)	Usable Length (cm)	Sheath Size (F)	Guidewire Acceptance (in)	Outside Diameter (mm)
SERP65-9-20-80	9	20	7.5-8.5	80	6	0.035	0.079
SERP65-9-30-80	9	30	7.5-8.5	80	6	0.035	0.079
SERP65-9-40-80	9	40	7.5-8.5	80	6	0.035	0.079
SERP65-9-60-80	9	60	7.5-8.5	80	6	0.035	0.079
SERP65-9-80-80	9	80	7.5-8.5	80	6	0.035	0.079
SERP65-10-20-80	10	20	8.5-9.5	80	6	0.035	0.079
SERP65-10-30-80	10	30	8.5-9.5	80	6	0.035	0.079
SERP65-10-40-80	10	40	8.5-9.5	80	6	0.035	0.079
SERP65-10-60-80	10	60	8.5-9.5	80	6	0.035	0.079
SERP65-10-80-80	10	80	8.5-9.5	80	6	0.035	0.079
SERP65-12-20-80	12	20	9.5-11.0	80	6	0.035	0.079
SERP65-12-30-80	12	30	9.5-11.0	80	6	0.035	0.079
SERP65-12-40-80	12	40	9.5-11.0	80	6	0.035	0.079
SERP65-12-60-80	12	60	9.5-11.0	80	6	0.035	0.079
SERP65-12-80-80	12	80	9.5-11.0	80	6	0.035	0.079
SERP65-9-20-120	9	20	7.5-8.5	120	6	0.035	0.079
SERP65-9-30-120	9	30	7.5-8.5	120	6	0.035	0.079
SERP65-9-40-120	9	40	7.5-8.5	120	6	0.035	0.079
SERP65-9-60-120	9	60	7.5-8.5	120	6	0.035	0.079
SERP65-9-80-120	9	80	7.5-8.5	120	6	0.035	0.079
SERP65-10-20-120	10	20	8.5-9.5	120	6	0.035	0.079
SERP65-10-30-120	10	30	8.5-9.5	120	6	0.035	0.079
SERP65-10-40-120	10	40	8.5-9.5	120	6	0.035	0.079
SERP65-10-60-120	10	60	8.5-9.5	120	6	0.035	0.079
SERP65-10-80-120	10	80	8.5-9.5	120	6	0.035	0.079
SERP65-12-20-120	12	20	9.5-11.0	120	6	0.035	0.079
SERP65-12-30-120	12	30	9.5-11.0	120	6	0.035	0.079
SERP65-12-40-120	12	40	9.5-11.0	120	6	0.035	0.079
SERP65-12-60-120	12	60	9.5-11.0	120	6	0.035	0.079
SERP65-12-80-120	12	80	9.5-11.0	120	6	0.035	0.079
SERP65-14-20-80	14	20	11.5-13.0	80	6	0.035	0.079
SERP65-14-30-80	14	30	11.5-13.0	80	6	0.035	0.079
SERP65-14-40-80	14	40	11.5-13.0	80	6	0.035	0.079
SERP65-14-60-80	14	60	11.5-13.0	80	6	0.035	0.079
SERP65-14-80-80	14	80	11.5-13.0	80	6	0.035	0.079
SERP65-14-20-120	14	20	11.5-13.0	120	6	0.035	0.079
SERP65-14-30-120	14	30	11.5-13.0	120	6	0.035	0.079
SERP65-14-40-120	14	40	11.5-13.0	120	6	0.035	0.079
SERP65-14-60-120	14	60	11.5-13.0	120	6	0.035	0.079
SERP65-14-80-120	14	80	11.5-13.0	120	6	0.035	0.079

Medtronic reserves the right to modify specifications without prior notice.

# Protégé™ RX

## Self-expanding Peripheral Stent System

Predictable deployment, visible results.

Reference Number	Unconstrained Stent Diameter (mm)	Unconstrained Stent Length (mm)	Rec. Lumen Size (mm)	Usable Catheter Length (cm)	Sheath Size (F)	Guidewire Acceptance
------------------	-----------------------------------	---------------------------------	----------------------	-----------------------------	-----------------	----------------------

### Straight

SEPX-6-20-135	6	20	4.5-5.5	135	6	0.014
SEPX-6-30-135	6	30	4.5-5.5	135	6	0.014
SEPX-6-40-135	6	40	4.5-5.5	135	6	0.014
SEPX-6-60-135	6	60	4.5-5.5	135	6	0.014
SEPX-7-20-135	7	20	5.5-6.5	135	6	0.014
SEPX-7-30-135	7	30	5.5-6.5	135	6	0.014
SEPX-7-40-135	7	40	5.5-6.5	135	6	0.014
SEPX-7-60-135	7	60	5.5-6.5	135	6	0.014
SEPX-8-20-135	8	20	6.5-7.5	135	6	0.014
SEPX-8-30-135	8	30	6.5-7.5	135	6	0.014
SEPX-8-40-135	8	40	6.5-7.5	135	6	0.014
SEPX-8-60-135	8	60	6.5-7.5	135	6	0.014
SEPX-9-20-135	9	20	7.5-8.5	135	6	0.014
SEPX-9-30-135	9	30	7.5-8.5	135	6	0.014
SEPX-9-40-135	9	40	7.5-8.5	135	6	0.014
SEPX-9-60-135	9	60	7.5-8.5	135	6	0.014
SEPX-10-20-135	10	20	8.5-9.5	135	6	0.014
SEPX-10-30-135	10	30	8.5-9.5	135	6	0.014
SEPX-10-40-135	10	40	8.5-9.5	135	6	0.014
SEPX-10-60-135	10	60	8.5-9.5	135	6	0.014

### Tapered

SEPX-8-6-30-135	8/6	30	(6.5-7.5)-(4.5-5.5)	135	6	0.014
SEPX-8-6-40-135	8/6	40	(6.5-7.5)-(4.5-5.5)	135	6	0.014
SEPX-10-7-30-135	10/7	30	(8.5-9.5)-(5.5-6.5)	135	6	0.014
SEPX-10-7-40-135	10/7	40	(8.5-9.5)-(5.5-6.5)	135	6	0.014

# Visi-Pro™

## Balloon-expandable Peripheral Stent System

Visible precision.

CATHETER LENGTH		STENT SIZE		BALLOON	COMPATIBILITY	
Reference Number 80 cm	Reference Number 135 cm	Diameter (mm)	Length (mm)	Balloon Length (mm)	Rec. Introducer Sheath (F)	Guidewire (in)
PXP35-05-12-080	-	5	12	15	6*	0.035"
PXP35-05-17-080	PXP35-05-17-135	5	17	20	6*	0.035"
PXP35-05-27-080	PXP35-05-27-135	5	27	30	6*	0.035"
PXP35-05-37-080	PXP35-05-37-135	5	37	40	6*	0.035"
PXP35-05-57-080	PXP35-05-57-135	5	57	60	6*	0.035"
PXP35-06-12-080	NA	6	12	15	6*	0.035"
PXP35-06-17-080	PXP35-06-17-135	6	17	20	6*	0.035"
PXP35-06-27-080	PXP35-06-27-135	6	27	30	6*	0.035"
PXP35-06-37-080	PXP35-06-37-135	6	37	40	6*	0.035"
PXP35-06-57-080	PXP35-06-57-135	6	57	60	6*	0.035"
PXP35-07-12-080	-	7	12	15	6*	0.035"
PXP35-07-17-080	PXP35-07-17-135	7	17	20	6*	0.035"
PXP35-07-27-080	PXP35-07-27-135	7	27	30	6*	0.035"
PXP35-07-37-080	PXP35-07-37-135	7	37	40	6*	0.035"
PXP35-07-57-080	PXP35-07-57-135	7	57	60	6*	0.035"
PXP35-08-17-080	PXP35-08-17-135	8	17	20	6*	0.035"
PXP35-08-27-080	PXP35-08-27-135	8	27	30	6*	0.035"
PXP35-08-37-080	PXP35-08-37-135	8	37	40	6*	0.035"
PXP35-08-57-080	PXP35-08-57-135	8	57	60	6*	0.035"
PXP35-09-17-080	PXP35-09-17-135	9	17	20	7	0.035"
PXP35-09-27-080	PXP35-09-27-135	9	27	30	7	0.035"
PXP35-09-37-080	PXP35-09-37-135	9	37	40	7	0.035"
PXP35-09-57-080	PXP35-09-57-135	9	57	60	7	0.035"
PXP35-10-17-080	PXP35-10-17-135	10	17	20	7	0.035"
PXP35-10-27-080	PXP35-10-27-135	10	27	30	7	0.035"
PXP35-10-37-080	PXP35-10-37-135	10	37	40	7	0.035"
PXP35-10-57-080	PXP35-10-57-135	10	57	60	7	0.035"

\*6 F = 0.085" I.D.